


TempTracer BLE - Start

1


- Login to novartis.ccis.se
- Select "Cloud Launch"



2

In Cloud Launch

- Select Profile (or scan barcode with profile information)
- Enter TripTags



3

Scan Logger ID(s) to be associated to shipment



4

Click "Configure"

ENTER LOGGER ID(S)

B1E2020110300150 x

Configure

5

Press and hold START button for 3 seconds, a green light will be shown, indicating that the logger has started



6

Press START button once if green light is shown the logger is started. If red, the logger is stopped

